

TEA-GSE DOCUMENTATION AND PROCEDURE

FRIDAY 18th, DECEMBER, 2017

MICAH JOHNSON

Introduction

This outlines the procedure for using the TEA GSE (Transfer ESIS Alignment Ground Station Equipment), as well as contains all the electrical and mechanical parts, custom part drawings, electrical schematics, codes, troubleshooting, and potential improvements. The instrument should be capable of position measurements of the ESIS diffraction gratings to a precision of 14 micron in three channels.

When run, TEA will move towards and away from a grating, measuring its position. The positions will be saved in a test report to be compared to when the alignment is being transferred.

Procedure

Startup

1. Turn all lasers and computer on
2. Once the computer is warmed up, open Pycharm.
3. Plug the USB for the Agilent into the computer, and turn it on.
4. Plug the USB for the K-cube into the computer, and turn it on.
5. Turn the voltage source on, and set it to 5 Volts.
6. Place thermal isolating box (like a thick box of tin-foil) over the stages and gratings, and tape it down to thermally isolate the system.
7. Wait approximately 3 hours for the lasers to warm up and for everything to come to thermal equilibrium.

Running

1. If the photodiodes are going to be needed, and you are not just doing a thermal measurement, turn the switches on the electronics box to on.
2. Edit the configurations for `vis_euv_trans.py` in Pycharm to [start position (mm), stop position (mm), number of positions to measure, number of times to loop through the positions, grating type (VIS or EUV), grating number, comment]. (Note the comment will not be saved).
3. Run the program. (This will take several hours).
4. Allow it to run until a pdf test report is created.
5. It may prompt you to enter a new configuration. If so, input a new start position, end position, number of steps, and number of measurements per step.

Parts

The following parts are all the parts used in the buildup of the TEA GSE that were not custom designed. Custom parts can be found the the Mechanical section of this document.

| Common Name | Part Number |
|--------------------|--------------------------------|
| Voltmeter | Agilent 34970A |
| Voltage Source | Agilent E3631A |
| Battery | Any 9 Volt battery |
| Laser | Thorlabs CPS635S |
| Laser Power Supply | Thorlabs LDS5 |
| Lens | Thorlabs LA1576-A |
| Beam splitter | Thorlabs BS007 |
| Retaining Right | Thorlabs SM9RR |
| Photodiode | Digikey Everlight 1080-1141-ND |
| Motor | Thorlabs MT1-Z8 |
| Motor Base Plate | Thorlabs MT401 |
| Motor Driver | Thorlabs KPS101 |
| Manuals Stage | Thorlabs PT1 |

Software

Overview

The software for finding the position of the gratings runs Python and Matlab. Matlab can be opened by directing to the path `/usr/local/MATLAB/R2017a/bin`, and running `./matlab`. Pycharm can be opened by directing to the `/Downloads/pycharm-community-2017.1.3/bin`, and running `./pycharm.sh`. Matlab does not need to be open to run the software, but is recommended if small changes are needed.

When the software runs, Python (`vis_euv_trans.py`) collects voltage measurements from the Agilent and moves the stepper motor. The voltage measurements from the thermistors are converted into temperature using the SteinhartHart equation (`temp_conversion.py`). The thermistor temperatures are then sent to Matlab where they will be plotted (`tmp.m`). Voltages from the photodiode and positions of the stepper motor are recorded and sent to Matlab. Matlab (`TEA_GSE_Final_Align.m`) will plot the voltages, output new start/stop/stepsizes for the motor, and test report.

There is also software from taking a series of temperature measurements without also driving the stepper motor (`long_temp.py`). It also uses the SteinhartHart equation (`temp_conversion.py`). Matlab then plots the temperature with time (`temp_long.m`). The configuration for the program is only the number of seconds that you want the temperature measurements taken over.

All software can be pulled from Gitlab at https://titan.ssel.montana.edu/gitlab/MOSES-ESIS/ESIS/Optics/Transfer_ESIS_Alignment_GSE.git. It also documents the installation procedure.

Matlab

TEA_GSE_Final_Align

```
1  
2  
3 %OVERVIEW: This program takes the data from the agilent  
for three  
4 %photodiode channels and outputs an adaptive grid which  
the stepper motor  
5 %and agelant will test over to get the precision of the  
peak intensity to
```

```

6 %be within a certain tolerance. A test report will be
7     generated after
8 %
9 %
10 % the following tripple commented lines are to be used
11 % for testing purposes
12 % by creating a simulated data matrix
13 % % % clear all
14 % % % close all
15 % % %
16 % % % stpsz=.002;
17 % % % iteration=1;
18 % % % dth=2.5:stpsz :5.5;
19 % % % vth1=-(dth-3.2).^2+dth+5;
20 % % % vth2=-(dth-3.1).^2+dth+4.5;
21 % % % vth3=-(dth-3.35).^2+dth+4.3;
22 % % %
23 % % % looops=10;
24 % % % n=0;
25 % % % dthex=[];
26 % % % vthex1=[];
27 % % % vthex2=[];
28 % % % vthex3=[];
29 % % %
30 % % % while n<looops
31 % % %     n=n+1;
32 % % %     dthex=[dthex dth];
33 % % %     vthex1=[vthex1 vth1];
34 % % %     vthex2=[vthex2 vth2];
35 % % %     vthex3=[vthex3 vth3];
36 % % % end
37 % % %
38 % % % dnos=((rand(1,(length(dth)*looops))-5)/5000)+1;
39 % % % vnos1=((rand(1,(length(dth)*looops))-5)/20)+1;
40 % % % vnos2=((rand(1,(length(dth)*looops))-5)/20)+1;
41 % % % vnos3=((rand(1,(length(dth)*looops))-5)/20)+1;
42 % % % vthnos1=vnos1.*vthex1;
43 % % % vthnos2=vnos2.*vthex2;
44 % % % vthnos3=vnos3.*vthex3;

```

```

45 % % % dthnos=dnos.*dthex;
46 %
47 % % % data=[dthnos;vthnos1;vthnos2;vthnos3]';
48 % % % csv_dir='/home/krg/Desktop/tester/iteration_3/
      data.csv';
49 % % % grating=3;
50 % % % type='VIS';

51
52
53
54 %input and outputs arguments
55 function [cont_loop, strt, stp, stept, mest, iteration]
    = TEA_GSE_Final_Align(csv_dir,stpsz,iteration,
                           grating,type)

56
57 %polynomial fits are normally not ideal
58 warning('off','all');

59
60 %data will be interpreted from its csv format
61 iteration=iteration+1;
62 data=csvread(csv_dir);
63 iteration=iteration+1;

64
65 %stage location at measurement
66 d=data(:,1);
67 v1=data(:,2);
68 v2=data(:,3);
69 v3=data(:,4);

70
71 strt=min(d);
72 stp=max(d);
73 stps=round((stpstrt)/(stpsz))+1;
74 lps=round(length(d)/(stps-1));

75
76 %loop prep
77 zerlps=zeros(stps,lps);
78 ds=zerlps;
79 vs1=zerlps;
80 vs2=zerlps;
81 vs3=zerlps;
82 %takes the distance vs voltage data array and creates a

```

```

        matrix of simaler
83 %distance measurements to be averaged together later ,
     and matracies of
84 %voltage measurements made at identical position to be
     bootstrapped later
85 for i=1:lps
86     ds (: , i)=d((( i-1)*stps+1):( i*stps));
87     vs1 (: , i)=v1((( i-1)*stps+1):( i*stps));
88     vs2 (: , i)=v2((( i-1)*stps+1):( i*stps));
89     vs3 (: , i)=v3((( i-1)*stps+1):( i*stps));
90 end
91
92 [ h , ^ ]= size ( vs1 );
93
94 %loop prep
95 zerstp=zeros(1,stps);
96 dav=zerstp;
97 vav1=zerstp;
98 vav2=zerstp;
99 vav3=zerstp;
100 %finds the average voltage for a given distance
101 for n=1:stps
102     dav(n)=mean( ds (n,:));
103     vav1(n)=mean( vs1(n,:));
104     vav2(n)=mean( vs2(n,:));
105     vav3(n)=mean( vs3(n,:));
106 end
107
108 %does a 4th order polynomial fit to all data
109 pd=min(d) : .0001 : max(d);
110
111 p1=polyfit(d,v1,4);
112 p2=polyfit(d,v2,4);
113 p3=polyfit(d,v3,4);
114
115 pruf1=polyval(p1,pd);
116 pruf2=polyval(p2,pd);
117 pruf3=polyval(p3,pd);
118
119 %finds over which data range to base fit
120 off1=1.5;

```

```

121 off2=1.75;
122 off3=2;
123
124 %does a polynomial fit over relavent data. The higher
      the iteration number,
125 %the wider the range over which is fit
126 if iteration==1
127     pd1=(dav(vav1==max(vav1))-off1):.0001:(dav(vav1==
           max(vav1))+off1);
128     pd2=(dav(vav2==max(vav2))-off1):.0001:(dav(vav2==
           max(vav2))+off1);
129     pd3=(dav(vav3==max(vav3))-off1):.0001:(dav(vav3==
           max(vav3))+off1);
130 end
131 if iteration==2
132     pd1=(dav(vav1==max(vav1))-off2):.0001:(dav(vav1==
           max(vav1))+off2);
133     pd2=(dav(vav2==max(vav2))-off2):.0001:(dav(vav2==
           max(vav2))+off2);
134     pd3=(dav(vav3==max(vav3))-off2):.0001:(dav(vav3==
           max(vav3))+off2);
135 end
136 if iteration==3
137     pd1=(dav(vav1==max(vav1))-off3):.0001:(dav(vav1==
           max(vav1))+off3);
138     pd2=(dav(vav2==max(vav2))-off3):.0001:(dav(vav2==
           max(vav2))+off3);
139     pd3=(dav(vav3==max(vav3))-off3):.0001:(dav(vav3==
           max(vav3))+off3);
140 end
141
142 dre11=dav(dav>=min(pd1) & dav<=max(pd1));
143 dre12=dav(dav>=min(pd2) & dav<=max(pd2));
144 dre13=dav(dav>=min(pd3) & dav<=max(pd3));
145 vrel11=vav1(dav>=min(pd1) & dav<=max(pd1));
146 vrel12=vav2(dav>=min(pd2) & dav<=max(pd2));
147 vrel13=vav3(dav>=min(pd3) & dav<=max(pd3));
148
149 pav1=polyfit(dre11,vrel11,4);
150 pav2=polyfit(dre12,vrel12,4);
151 pav3=polyfit(dre13,vrel13,4);

```

```

152
153 pavd1=polyval(pav1, pd1);
154 pavd2=polyval(pav2, pd2);
155 pavd3=polyval(pav3, pd3);

156
157 %plots all data with fits
158 all_fig_combo=figure(1);
159 hold on
160 p1=scatter(d,v1, 'r', 'MarkerFaceAlpha', .5, ...
    'MarkerEdgeAlpha', 0.05);
161 p2=scatter(d,v2, 'b', 'MarkerFaceAlpha', .5, ...
    'MarkerEdgeAlpha', 0.05);
162 p3=scatter(d,v3, 'g', 'MarkerFaceAlpha', .5, ...
    'MarkerEdgeAlpha', 0.05);
163 p4=scatter(dav,vav1, 'm', 'MarkerFaceAlpha', .5, ...
    'MarkerEdgeAlpha', 0.1);
164 p5=scatter(dav,vav2, 'c', 'MarkerFaceAlpha', .5, ...
    'MarkerEdgeAlpha', 0.1);
165 p6=scatter(dav,vav3, 'k', 'MarkerFaceAlpha', .5, ...
    'MarkerEdgeAlpha', 0.1);
166 p7=plot(pd, pruf1, 'r');
167 p8=plot(pd, pruf2, 'b');
168 p9=plot(pd, pruf3, 'g');
169 p10=plot(pd1, pavd1, 'm');
170 p11=plot(pd2, pavd2, 'c');
171 p12=plot(pd3, pavd3, 'k');
172 legend([p7 p8 p9 p10 p11 p12], {'All Data 1', 'All Data 2', ...
    'All Data 3', 'Averaged Data 1', 'Averaged Data 2', ...
    'Averaged Data 3'})
173 xlabel('Distance (mm)')
174 ylabel('Voltage (V)')
175 title('Data and Averaged Data with Fits')

176
177 %saves as png. pdf would be preferable, but file is too
    large it crashes
178 %computer
179 print(all_fig_combo, strrep(csv_dir, 'data.csv', ...
    all_fig_combo), '-dpng')

180
181 %finds peak of data
182 dmvl=dav(vav1==max(vav1));

```

```

183 dmv2=dav(vav2==max(vav2));
184 dmv3=dav(vav3==max(vav3));
185
186 %adaptive grid for 2nd run through
187 if iteration==1
188     strt=min([dmv1 dmv2 dmv3])-off1;
189     stp=max([dmv1 dmv2 dmv3])+off1;
190     mest=20;
191     stept=.002;
192     cont_loop=0;
193 end
194
195 %adaptive grid for 3rd run through. stops if 4th run
196 %through. outputs pdf
197 %test report.
198 if iteration==2 || iteration==3
199     %loop prep
200     times=1000;
201     tim=0;
202     mxxs1=zeros(1,times);
203     mxxs2=zeros(1,times);
204     mxxs3=zeros(1,times);
205     vbot1=zeros(1,lps);
206     vbot2=zeros(1,lps);
207     vbot3=zeros(1,lps);
208
209     while(tim<times)
210
211         %display progress through loop
212         if rem(10*tim,times)==0
213             sprintf('>%g%%',round(100*tim/times))
214         end
215         tim=tim+1;
216
217         %samples with replacement from measurements at
218         %every distance
219         for i=1:h
220             vbot1(i)=mean(datasample(vs1(i,:),lps));
221             vbot2(i)=mean(datasample(vs2(i,:),lps));
222             vbot3(i)=mean(datasample(vs3(i,:),lps));
223         end

```

```

222
223 %does polynomaial fit to bootstrapped data
224 pbot1=polyfit(dav,vbot1,4);
225 pbot2=polyfit(dav,vbot2,4);
226 pbot3=polyfit(dav,vbot3,4);
227 ft1=polyval(pbot1,pd1);
228 ft2=polyval(pbot2,pd2);
229 ft3=polyval(pbot3,pd3);
230 mx1=pd1(ft1==max(ft1));
231 mx2=pd2(ft2==max(ft2));
232 mx3=pd1(ft3==max(ft3));
233 %find max for bootstrapped data
234 mxxs1(tim)=mx1(1);
235 mxxs2(tim)=mx2(1);
236 mxxs3(tim)=mx3(1);
237 end
238
239 %creates a confidence interval
240 mxxs1=sort(mxxs1);
241 mxxs2=sort(mxxs2);
242 mxxs3=sort(mxxs3);
243
244 num_pt=round(times*.98/2,0);
245
246 lwd1=mxxs1(times/2-num_pt);
247 lwd2=mxxs2(times/2-num_pt);
248 lwd3=mxxs3(times/2-num_pt);
249 hbd1=mxxs1(times/2+num_pt);
250 hbd2=mxxs2(times/2+num_pt);
251 hbd3=mxxs3(times/2+num_pt);
252
253 %plots histograms with 98% confidence intervals
254 hist1=figure();
255 hold on
256 hist(mxxs1,times/10)
257 line([lwd1 lwd1],ylim,'Color','r');
258 line([hbd1 hbd1],ylim,'Color','r');
259 hist2=figure();
260 hold on
261 hist(mxxs2,times/10)
262 line([lwd2 lwd2],ylim,'Color','r');

```

```

263     line([ hibd2 hibd2] , ylim , 'Color' , 'r') ;
264     hist3=figure();
265     hold on
266     hist(mxxs3 , times/10)
267     line([ lwd3 lwd3] , ylim , 'Color' , 'r') ;
268     line([ hibd3 hibd3] , ylim , 'Color' , 'r') ;
269
270 %convert to micron
271 rng1=(hibd1-lwd1)*1000;
272 rng2=(hibd2-lwd2)*1000;
273 rng3=(hibd3-lwd3)*1000;
274
275 %largest confidence interval
276 hi_rng=max([rng1 rng2 rng3]);
277
278 %change according to resolution buget. adaptive
279 %grid for 3rd run
280 %through.
281 if hi_rng<=15.3
282     pf='Passed';
283     strt=min([dmv1 dmv2 dmv3])-off2;
284     stp=max([dmv1 dmv2 dmv3])+off2;
285     mest=30;
286     stept=.002;
287     cont_loop=0;
288 else
289     pf='Failed';
290     strt=0;
291     stp=0;
292     mest=0;
293     stept=0;
294     cont_loop=2;
295 end
296
297 %saves histograms. pdf would be better
298 print(hist1,strrep(csv_dir,'data.csv','Histogram 1','-dpng'))
299 print(hist2,strrep(csv_dir,'data.csv','Histogram 2','-dpng'))
300 print(hist3,strrep(csv_dir,'data.csv','Histogram 3','-dpng'))

```

```

300
301 %begin formatting for pdf output
302 %all fig combo file location
303 alfg=strrep(csv_dir,'data.csv','all_fig_combo.png')
304 ;
305 %variables from script to be used in test report
variableList={type num2str(grating) num2str(mode(
    mxxs1)) num2str(rng1) num2str(mode(mxxs2))
    num2str(rng2) num2str(mode(mxxs3)) num2str(rng3)
    num2str(iteration) num2str(lps) num2str(stps)
    alfg pf};

306
307 %tex lines with script inputs
308 line0 = cellstr(['{\textbf{Test}',variableList{13},
    '} \par']);
309 line1 = cellstr(['{\textsf{large Grating Type:}',%
    variableList{1}, '\par']]);
310 line2 = cellstr(['{\textsf{large Grating Number:}',%
    variableList{2}, '\par']]);
311 line3 = cellstr(['\indent... for channel 1 lies at',
    'variableList{3} ' mm to a precision of %
    variableList{4} ' microns.\n']);
312 line4 = cellstr(['\indent... for channel 2 lies at',
    'variableList{5} ' mm to a precision of %
    variableList{6} ' microns.\n']);
313 line5 = cellstr(['\indent... for channel 3 lies at',
    'variableList{7} ' mm to a precision of %
    variableList{8} ' microns.\n']);
314 line6 = cellstr(['\indent... after ' variableList{9}
    ' intensity scans.\n']);
315 line7 = cellstr(['\indent... with ' variableList{10}
    ' measurements taken at each distance.\n']);
316 line8 = cellstr(['\indent... with data from',
    'variableList{11} ' distances.\n']);
317 line9 = cellstr(['\includegraphics[]{' variableList
    {12} '}\n']);

318
319 %entire tex file
320 str = ['\documentclass[12pt,a4paper]{article}'
    '\usepackage{datetime}'
    '\usepackage{ragged2e}'
```

```

323   '\usepackage{graphicx}'  

324   '\usepackage{subcaption}'  

325   '\usepackage{mwe}'  

326   '\usepackage{float}'  

327   '\usepackage{grffile}'  

328   '\pagenumbering{gobble}'  

329   ,,  

330   '\begin{document}'  

331   ,,  

332   '\begin{center}'  

333   '{\scshape\LARGE VIS-EUV Test Report \par}'  

334   '{\scshape\Large \today, \currenttime \par}'  

335   '\bigskip'  

336   line1  

337   line2  

338   '\bigskip'  

339   line0  

340   '\end{center}'  

341   ,,  

342   '\noindent We are 98\% confident that the  

343   peak intensity ...\\'  

344   line3  

345   line4  

346   line5  

347   ,  

348   '\noindent The final model was created...\\'  

349   line6  

350   line7  

351   line8  

352   ,  

353   '\begin{figure}[H]'  

354   '\centering'  

355   line9  

356   '\end{figure}'  

357   ,  

358   '\end{document}'];  

359  

360 %file path of tex document  

361 pth=csv_dir(1:end-20);  

362 pthfol=strcat(pth,'reports');  

363 pth=strcat(pth,'reports/Test_Report.tex');

```

```

363
364 %prepare file for tex report
365 formatSpec = '%s\r\n'; %specify the data type
366 % being output through fprintf
367 [nRows,~] = size(str); %tell matlab the size of
368 % the file
369 mkdir(pthfol);
370 fileID = fopen(pth, 'wt'); %open the blank report
371 %file
372
373 %writes tex file line by line
374 for row = 1:nRows
375     fprintf(fileID,formatSpec,str{row,:}); %%
376         print the data onto the file
377 end
378
379 fclose(fileID); %close the finished file
380 cd (pth(1:end-15));
381
382 %compile and open latex pdf document
383 system('pdflatex Test_Report.tex > /dev/null 2>&1')
384 ;
385 system('xdg-open Test_Report.pdf');
386
387 end
388
tmp

1
2
3 function tmp(csv_dir)
4
5 dat=csvread(csv_dir);
6
7 sec=dat(:,16);
8
9 temp=figure();
10 hold on
11 for i=[1:4 6:10]
12     t=dat(:,(i+5));
13     plot(sec,t)

```

```

14 end

15
16 legend( '1' , '2' , '3' , '4' , '6' , '7' , '8' , '9' , '10' )
17 xlabel( 'Time ( seconds )' )
18 ylabel( 'Temperature (F)' )
19 title( 'Temperature in Time' )

20
21 print( temp , strrep( csv_dir , 'data.csv' , 'Temperature' ) , '-
dpdf' );

temp_long

1
2
3 function temp_long( csv_dir )
4
5 dat=csvread( csv_dir );
6
7 sec=dat (:,12);
8
9 temp=figure();
10 hold on
11 for i=[1:4 6:10]
12     t=dat (:,( i+1));
13     plot(sec ,t)
14 end
15
16 legend( '1' , '2' , '3' , '4' , '6' , '7' , '8' , '9' , '10' )
17 xlabel( 'Time ( seconds )' )
18 ylabel( 'Temperature (F)' )
19 title( 'Temperature in Time' )

20
21 print( temp , strrep( csv_dir , 'data.csv' , 'Temperature' ) , '-
dpdf' );

```

Python

vis_euv_trans.py

```

1 # Python 3 Program for moving ThorLabs Z812B motorized
   stage and taking intensity measurements with
   Agilent 34970A

```

```

2 # Calling sequence: python3 vis_euv_trans.py [start
   position] [end position] [number of steps] [number
   of measurements per step]
3
4 import thorpy
5 import thorpy.message
6 import thorpy.comm.discovery as comm
7 import thorpy.stages
8 import agilent
9 from temp_conversion import steinhart
10 import os
11 import csv
12 import time
13 from time import sleep
14 import sys
15 import matlab.engine
16 eng = matlab.engine.start_matlab()
17 eng.cd(r'../Matlab')
18
19 # Check if the program received the correct number of
   arguments
20 if len(sys.argv) != 8:
21     sys.exit("Incorrect number of arguments! \n Calling
       sequence: python3 vis_euv_trans.py [start
       position] [end position] [number of steps] [
       number of measurements per step] [VIS | EUV] [
       Grating #] [Comment] \n" + str(sys.argv))
22
23 # Set up measurement variables
24 start_pos = float(sys.argv[1])      # Position where the
   sequence starts, measured from the home position
25 end_pos = float(sys.argv[2])        # Position where the
   sequence ends, measured from the home position
26 num_steps = int(sys.argv[3])        # Number of steps to
   take with the motor between start and end position
27 num_meas = int(sys.argv[4])         # Number of intensity
   measurements to take per step
28
29 # Type of grating to analyze for this measurement (
   Visible or EUV)
30 grating_type = str(sys.argv[5])

```

```

31 if(grating_type != "VIS" and grating_type != "EUV"):
32     sys.exit("Incorrect grating type argument, use 'VIS'
33             ' or 'EUV'")
34 # Serial number of the grating under test
35 grating_num = int(sys.argv[6])
36 if(grating_num < 1 or grating_num > 6):
37     sys.exit("Incorrect grating number specified. Use a
38             value between 1 and 6")
39 # Comment associated with this test run
40 grating_comment = str(sys.argv[7])
41
42 # prepare the directory structure
43 grating_dir = "/home/krg/ESIS_VisEUV_transfer_sw/Output
44             /" + grating_type + "/" + "grating_" + str(
45             grating_num)
46 if not os.path.exists(grating_dir):
47     os.makedirs(grating_dir)
48
49 timestr = time.strftime("%Y%m%d-%H%M%S")
50 exp_dir = grating_dir + "/" + timestr
51 os.makedirs(exp_dir)
52
53 report_dir = exp_dir + "/reports"
54 os.makedirs(report_dir)
55 iter_dir = exp_dir + "/iterations"
56 os.makedirs(iter_dir)
57
58 cont_loop = 0      # 0 = for continuing adptive mesh
59             refinement, 1 = for manual mesh refinement, 2 =
60             bootstrapping converged, exit program.
61 iteration = 0
62
63 # Open the Z812B motorized stage
64 motor = None
65 x = comm.discover_stages()
66 for i in x:
67     motor = i

```

```

66
67 if motor == None:
68     sys.exit("Z812 stage not connected")
69
70 # Open the Agilent 34970A for voltage measurements
71 voltmeter = agilent.init_serial()
72
73 # Prepare the motor for a sequence
74 motor.acceleration = 0.5
75 motor.max_velocity = 2.3      # Set max velocity
76                         parameter on Z812 (mm/s)
77 motor.print_state()
78
79 motor.home()
80
81 while(True):
82
83     # Check the return value of Matlab script to
84     # determine what to do next
85     if(cont_loop == 1):
86         [start_pos, end_pos, num_steps, num_meas] =
87             input("Adaptive grid selection failed,
88                   Please enter the new gridsize... ")
89     elif(cont_loop == 2):
90         break
91
92     dx = (end_pos - start_pos) / num_steps    # Amount to
93     move for each step
94
95     # Check that the arguments are physically
96     # attainable
97     if start_pos == end_pos:
98         sys.exit("Nothing to do. Start position equal
99                 to end position")
100
101    if min(start_pos, end_pos) < motor.
102        home_offset_distance:
103            sys.exit("Minimum value too close to end of
104                  motor")
105
106    if max(start_pos, end_pos) > (12.0 - motor.

```

```

    home_offset_distance):
98     sys.exit("Maximum value too close to end of
      motor")
99
100    if num_meas <= 0:
101        sys.exit("Number of measurements must be larger
      than zero")
102
103    if num_steps <= 0:
104        sys.exit("Number of steps must be larger than
      zero")
105
106    # Prepare directory for this iteration
107    cur_dir = iter_dir + "/iteration_" + str(iteration)
108    os.makedirs(cur_dir)
109
110    # Open CSV file for writing output
111    csv_fn = cur_dir + "/data.csv"
112    csvfile = open(csv_fn, 'w', newline=',')
113    csv_writer = csv.writer(csvfile)
114
115
116    # Loop through the sequence
117    for num in range(0, num_steps+1):
118
119        # Make sure the motor isn't going to be damaged
120        # by the move
121        if(motor.status_motor_current_limit_reached or
           motor.status_motion_error or motor.
           status_forward_hardware_limit_switch_active
           or motor.
           status_reverse_hardware_limit_switch_active):
122            sys.exit("Z812B encountered an error")
123
124        # Move the motor to the new position and wait
125        # for completion
126        motor.position = start_pos + num * dx
           sleep(0.2)
           while (motor.status_in_motion_forward or motor.
           status_in_motion_reverse or motor.

```

```

status_in_motion_jogging_forward or motor.
status_in_motion_jogging_reverse or motor.
status_in_motion_homing):
    sleep(0.1)
    #motor.print_state()

# Take the intesity measurements
data = []
data.append(str(motor.position))
for j in range(0, num_meas):           # Take
    intensity_measurements
    raw = agilent.measure(voltmet, 101, 101)
    rawlist = raw.split(',')
    for k in range(0, len(rawlist)):
        if k < 3:      # photodiode measurement
            data.append(rawlist[k])
        else:
            data.append(steinhart(rawlist[k]))

csv_writer.writerow(data)    # write next row of
                            data to file
 csvfile.flush()

 csvfile.close() # Close CSV file

[cont_loop, start_pos, end_pos, num_steps, num_meas,
 , iteration] = eng.complete_align_single(csv_fn,
                                           iteration, grating_num, grating_type, nargout=6)
cont_loop = int(cont_loop)
num_steps = int(num_steps)
num_meas = int(num_meas)
iteration= int(iteration)
print(cont_loop, start_pos, end_pos, num_steps,
      num_meas, iteration)

input("Press Enter to quit...")
```

```

temp_conversion

1 import math
2
3 def steinhart(V):
4     Rb = 4990
5     C1 = 1.285e-3
6     C2 = 2.362e-4
7     C3 = 9.285e-8
8     Vin = 12.
9
10 if V <= 0:
11     V = 1e-6
12
13 # R = -Rb * V / (V-Vin)
14 R = Rb * (Vin/V - 1)
15 return 1 / (C1 + C2 * math.log(R) + C3 * math.pow(
    math.log(R), 3)) - 273.15

long_temp

1
2
3 import thorpy
4 import thorpy.message
5 import thorpy.comm.discovery as comm
6 import thorpy.stages
7 import agilent
8 from temp_conversion import steinhart
9 import os
10 import csv
11 import time
12 from time import sleep
13 import sys
14 import matlab.engine
15 eng = matlab.engine.start_matlab()
16 eng.cd(r'../Matlab')
17
18 duration = float(sys.argv[1])
19
20
```

```

21 temp_dir = "/home/krg/Transfer_ESIS_Alignment_GSE/
22   control_software/Output/Temperature"
23 if not os.path.exists(temp_dir):
24     os.makedirs(temp_dir)
25
26 timestr = time.strftime("%Y%m%d-%H%M%S")
27 exp_dir = temp_dir + "/" + timestr
28 os.makedirs(exp_dir)
29
30 # Open the Agilent 34970A for voltage measurements
31 voltmet = agilent.init_serial()
32
33 csv_fn = exp_dir + "/data.csv"
34 csvfile = open(csv_fn, 'w', newline=',')
35 csv_writer = csv.writer(csvfile)
36
37 start_time=time.time()
38
39 dif = 0
40
41 while dif < duration:
42     data = []
43     raw = agilent.measure(voltmet,104,114)
44     rawlist = raw.split(',')
45     for k in range(0,(len(rawlist)+1)):
46         if k == 0: # photodiode & Vin measurement
47             data.append(rawlist[k])
48         elif k == 11:
49             dif = time.time()-start_time
50             data.append(dif)
51         else:
52             data.append(steinhart(float(rawlist[k]),
53                                   float(rawlist[0])))
54     csv_writer.writerow(data) # write next row of data
55       to file
56     csvfile.flush()
57 eng.temp_long(csv_fn, nargout=0)
58

```


⁵⁹ input("Press Enter to quit . . .")

Electrical


Overview

There are two electronic systems used in measurements for the TEA GSE. One is a simple voltage divider with a thermistor and resistor. This is used to measure the temperature. There are 10 repetitions of the temperature circuit so the temperature of 10 locations can be measured simultaneously. One of the channels does not work, and is labeled as such. The other circuit is meant for measuring the photodiodes. The photodiode circuit converts the current measurement of the photodiode into voltage, and filters out the high frequency noise. There are 3 repetitions of the photodiode circuit for each of the three channels. All outgoing wires are connected by D-subs.

Photodiode




Temperature





Mechanical


The parts designed and built for the grating alignment can be divided up into "TEA GSE", "External", and "Electrical". The TEA GSE includes the parts used to actually make the grating measurements. It can be attached to another kind of buildup to make measurements of another kind of mirror or grating, attached in any other way. The external buildup is used to hold the TEA, stages, and gratings. It is specific to the ESIS gratings held by the tuffet. The Electrical buildup holds all the electronics, apart from the voltmeter and voltage supply. The housing for the measurement electronics is also included.


TEA GSE












| REV. | | DESCRIPTION | | DATE | APPROVED |
|------|--|-----------------|--|-----------|------------|
| 1.0 | | INITIAL RELEASE | | 3/28/2017 | R. MEUCHEL |



D

C

B

A

REVISIONS

REV.
1.0
DESCRIPTION
INITIAL RELEASE
DATE
3/28/2017
APPROVED
R. MEUCHEL

TITLE:
Left_Right Laser Holder

REV
1.0
03-00024

SHEET 1 OF 1
SCALE: 2:1
DO NOT SCALE DRAWING
PROPRIETARY AND CONFIDENTIAL
DATE 3/28/2017
ENGINEER R. MEUCHEL


REV 1.0

UNLESS OTHERWISE SPECIFIED:
DIMENSIONS ARE IN INCHES
TOLERANCES:
ANGULAR +0.5°
XXX ± 0.020"
XXXXX ± 0.005"
XXXXXX ± 0.0005"

QUANTITY: 2

MATERIAL 6061-T6 (SS)
FINISH De-Burr
SHEET 1 OF 1
SCALE: 2:1
DO NOT SCALE DRAWING
PROPRIETARY AND CONFIDENTIAL
DATE 3/28/2017
ENGINEER R. MEUCHEL

1 2 3 4 5 6 7 8




D




C

B

A





| REV. | | DESCRIPTION | | DATE | APPROVED |
|------|--|-----------------|--|-----------|------------|
| 1.0 | | INITIAL RELEASE | | 3/28/2017 | R. MEUCHEL |





| | | | |
|--|--|--------------------------------------|------------------------|
| UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR $\pm 0.5^\circ$ $X_{XXX} \pm 0.020^\circ$ $X_{XXXX} \pm 0.005^\circ$ $X_{XXXXX} \pm 0.0005^\circ$ | | TITLE: Bottom Laser Holder | |
| 3RD ANGLE PROJECTION | | ENGINEER: R. MEUCHEL | REV. 1.0 |
| SIZE B | | PART NO. 03-00025 | SCALED SHEET 1 OF 1 |
| DATE 3/28/2017 | | DO NOT SCALE DRAWING | SCALE: 2:1 |
| PROPRIETARY AND CONFIDENTIAL | | | |

| MATERIAL | 6061-T6 (SS) | FINISH | De Burr |
|----------|--------------|--------|---------|
| 8 | 9 | 4 | 5 |












| REV. | | | | DESCRIPTION | | | | DATE | | APPROVED | |
|------|--|--|--|------------------|--|--|--|--------|--|------------|--|
| 1.0 | | | | Original Release | | | | 8/8/17 | | M. Johnson | |
| 1.1 | | | | | | | | | | | |
| 1.2 | | | | | | | | | | | |





| | | | | |
|------------------------------|---------------|---|------------|--------------|
| UNLESS OTHERWISE SPECIFIED: | | DIMENSIONS ARE IN INCHES | | |
| | | TOLERANCES: ANGULAR: $\pm 0.5^\circ$ $X.XXX \pm 0.020$ $X.XXXX \pm 0.005$ $X.XXXX \pm 0.0005$ | | |
| | | TITLE: Third_Apet | | |
| 3RD ANGLE PROJECTION | | | | |
| ENGINEER | M. Johnson | SIZE | PART NO. | REV |
| DATE | 8/8/17 | B | | 1.0 |
| PROPRIETARY AND CONFIDENTIAL | | DO NOT SCALE DRAWING | SCALE: 7.0 | SHEET 1 OF 1 |
| MATERIAL | Aluminum 5056 | | | |
| FINISH | Black Anodize | | | |


External Buildup














3.750

2.500

1.250

.500

.000



4.250

5.000

5.750

7.250

13X ϕ .089
UNC 4-40 TAP

.375

10.250

SCALE 0.500

Aluminum 5056



Original Release

8/10/17

M. Johnson

Tab_Adap

Deburr

8/10/17

M. Johnson

1.0

| REV. | | DESCRIPTION | | DATE | APPROVED |
|------|--|------------------|--|---------|------------|
| 1.0 | | Original Release | | 8/10/17 | M. Johnson |
| 1.1 | | | | | |
| 1.2 | | | | | |



4X ϕ .386 ∇ .26

.000
.500
1.500
2.500
3.500
4.500
5.000

4X ϕ .257 THRU

4X ϕ .120 THRU

REAM .1247 THRU

6X ϕ .189 THRU

UNC 1/4-20 THRU

8.000
8.500
9.500
10.000

11.500
12.000

SCALE 0.500

UNLESS OTHERWISE SPECIFIED:

DIMENSIONS ARE IN INCHES
TOLERANCES:
ANGULAR: $\pm 0.5^\circ$
 $XXX \pm 0.020$
 $XXXX \pm 0.005$
 $XXXXX \pm 0.0005$

TITLE: Tab_Adap

3RD ANGLE PROJECTION

ENGINEER: M. Johnson

DATE: 8/10/17

PART NO. B

REV. 1.0

SHEET 1 OF 1

MATERIAL Aluminum 5056

FINISH Deburr

SCALE: 1

2

3

4

5

6

7

8

D

C


B

D


C

B

A



Electrical





Examples




Figure 1: The temperature of each of the thermistor channels over 16 hours. The Lasers were turned on in the first 10 seconds of plotting. TEA was in a box of aluminum foil taped to the table throughout the duration of the test.

VIS-EUV TEST REPORT

WEDNESDAY 30TH AUGUST, 2017, 14:32

GRATING TYPE:VIS
GRATING NUMBER:1

We are 98% confident that the peak intensity...
...for channel 1 lies at 4.2401 mm to a precision of 952.5 microns with a target of 15.3 micron.
The final model was created...
...after 4 intensity scans.
...with 5 measurements taken at each distance.
...with data from 3465 distances.




Figure 2: An example test report for a single channel being read. A test report for all three channels will contain the same information.




Figure 3: This is the type of plot that should be created when a test is run. Work will continue to sharpen the intensity peaks.

Troubleshooting

1. Q — The programs are throwing unexpected errors, and the K-cube has stopped responding

A — Turn the Agilent and the K-cube off and unplug them from the computer. Make sure the computer is turned on, and logged in. First plug the Agilent USB in. Turn the Agilent on. Wait a few seconds and plug the K-cube into the computer. Turn the K-cube on. This should reset the way the computer reads from the USB's

2. Q — The power button of the K-cube has stopped working.

A — Unplug the K-cube, turn off the cube, plug it in, and turn it on. It tends to do then when the USB order has been messed up, so see the above answer.

3. Q — The program has stopped running after the first iteration

A — It is likely that one of the photodiode channel is not receiving enough light, and checking the graph should reveal that the curvature of one of channels is not very strong. This will mess with the adaptive grid algorithm. Provide more light to the channel in question, and the algorithm should be able to work properly.

Improvements

- Place an opto-isolator between the beam splitter and laser to prevent feedback into the laser to change intensity and laser modes.
- Connect an oscilloscope in parallel with a laser to check the stability of the laser in the system.
- Use a block machined to a greater precision than of the stepper motor to zero the motor before every test.
- Make all parts of titanium to prevent thermal fluctuations from influencing buildup.
- Test different diffusing glass between the beam splitter and photodiode to increase sensitivity.

- Replace amplifying resistor in photodiode circuit to increase gain and voltage output.
- Wire directly into Agilent card.
- Stop down aperture between the beam splitter and laser so the entire lens can be used.
- Create a system for organizing the thermistor wires.

Images

